



SAIKANG/Medical

Healthier your life

B4w

FEATURES AND HIGHLIGHTS











Technical parameters

Size (LxW)	2180x960 mm
 Bed height adjustment 	460-710 mm
 Mattress platform 	1925x830 mm
Safe working load	250 kg
Back-rest adjustment	0-75° (±10°)
Knee-rest adjustment	0-35° (±10°)
Reverse trendelenburg	0-12°
 Trendelenburg 	0-12°

Technical configuration

✓ Manual crank	4 pcs
☑ 5" hardshell castors with break	4 pcs
☑ Foldable side rail with ProErgo-Handle™	1 set
☑ ABS bed ends with safe lock	1 set
☑ IV pole	1 set
☑ IV pole prevision	4 pcs
☑ Drainage hook	2 pcs









Design

Use new and strong inserted ABS fireproof panel avoid off. Easy to be cleaned and installed.



Steels | Painting

Baosteel from Fortune Global 500. 11 Process epoxy painting, ASTM testing anti-baterial, paint thickness 0.12mm, brightness 60°, paint can resist 50kg impact.



Welding | Chassis

Panasonic robotic ensure, 360° full smooth wedlding Solid bed base made of 5.0 mm stamped rotating plate.



Hand Crank

Retractable hand crank, 100,000 times remaining intact.
Extreme protection device, sounds when tipped to the max position.



Wheels

TPR tire no worn out after running 30KM, save anti-winding hard shell, united forming without bolts.
Pass Dynamic test: bearing 120kg run 30KM and pass obstacles 500 times



Head & Foot Board

Full ABS material, head and foot board can be reversable.



Bed Surface

Stamping forming platform, smooth surface no welding scar, bearing 250KG load. Option ABS platform and resin X-ray translucent platform.



Side Rail

L-shaped foldable guardrail(Patent), bended tube, solid installing, 0 shaking, anti-pinch, ensure patient's safety.



IV Pole

Stainless steel solumn with plastic hook.



Drainage Hooks

Mobile plastic drainage hooks.

















Design

Research

Manufacturing

Maintenance

Service



USA

Saikang Medical Technology, Inc. 17800 Castleton St Ste 180, City of Industry, 91748, California. Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC Truda 174, 454080 Chelyabinsk. Tel.: +7 (351) 2145566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd. Viraj Complex Unit 3, Mombasa Road, P.O. Box 14262-00800, Nairobi. Kenya. Tel.: +254733704690

Middle East

Saikang Medical Supplies LLC X22-S03 England Cluster, International City, Dubai. United Arab Emirates. Tel.: +20 128 8962497









www.saikangmedical.com







Jiangsu Saikang Medical Equipment Co., Ltd. STOCK CODE: 870098

Company video

General catalogue